

Product Code . ILS-PLE-10131

Ring & Ball Apparatus






Description

Ring & Ball Apparatus

The ball expands due to increased molecular motion. When the Ring and Ball both are at same temperature the ball easily passes through the ring. Indian Lab Suppliers leading a Ring & Ball Apparatus is used for demonstrating thermal expansion of metal. Heat the ball in a flame. Try to pass the ball through the ring, it will not pass. The Ring & Ball Apparatus is a physics lab equipment used to teach heat, temperature and thermodynamics.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com