


## Indian Lab Suppliers









## **Bell Experiment Apparatus**

## **Description**

## **Bell Experiment Apparatus**

It consists of a 2 litre bell jar with a rubber stopper fitted with binding posts and connectors from which an electric bell is suspended. It operates on 3 to 6 volts, using either battery or a transformer. For demonstrating that sound cannot be transmitted in a vacuum.

Indian Lab Suppliers,

Direct Contact Details ← +91-8569909696 ☑ sales@indianlabsuppliers.com ☐ www.indianlabsuppliers.com