



Product Code . ILS-BLE-10367

Potometer Apparatus White Stand

Description

Potometer Apparatus White Stand

The apparatus can also be used to measure humidity by connecting a porous pot.

Used to show transpiration which occurs in small plant shoots. The apparatus is submerged and the specimen placed in the tube.




When transpiration occurs, the bubble within the capillary moves. A syringe is used to zero the scale.

Allows the students to measure transpiration with minimal risk of breakages.

The apparatus has only one glass component (the capillary tube) and is mounted on a strong backing plate.

Dimensions: 195 x 150 x 75mm.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com