

Product Code . ILS-BLE-10368

Ganong Potometer Apparatus



Description

Ganong Potometer Apparatus

The open end of the capillary tube is dipped into a beaker of water.

For demonstration of transpiration rates and rates of water absorption by a cut shoot under various external conditions.

A plant shoot is sealed into a stopper which fits the end tube and the whole apparatus is filled with water.

The position of the air bubble may easily be adjusted by opening the glass stopcock to admit water from the reservoir.

The whole apparatus is mounted on a varnished wooden base by spring clips.

Base: 20 x 80mm.

Capillary Bore: 1mm.

Contact jLab for your Educational School Science Lab Equipments. We are best school science lab equipments suppliers in india, school science lab instruments, school science lab instruments exporter, school science lab instruments manufacturer, school scientific equipments, science educational equipment in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com