



**Product Code . ILS-EAST-10850**

## Viscosity Apparatus Stokes Method

### Description

---

On the lower end a rubber cork without holes is fitted whereas on upper end a rubber cork with a funnel is provided.

With two rubber band index mark and 6 steel balls of different diameter.




Without lead shot.

A Borosilicate glass tube of about 30 cm diameter and one meter long is made to stand vertically firmly in between two upright bars with the help of clamps on a heavy cast iron base.

Both sides of the tube are open.

---

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  [sales@indianlabsuppliers.com](mailto:sales@indianlabsuppliers.com)  
 [www.indianlabsuppliers.com](http://www.indianlabsuppliers.com)