



Product Code . ILS-POM-10867

Bell Experiment Apparatus

Description

It consists of a 2 litre bell jar with a rubber stopper fitted with binding posts and connectors from which an electric bell is suspended.




It operates on 3 to 6 volts, using either battery or a transformer.

For demonstrating that sound cannot be transmitted in a vacuum.

To be used with Air Pump (Vacuum Pump with Aluminium Plate).

Contact jLab for your Educational School Science Lab Equipments. We are best technical educational equipments manufacturers, technical educational equipments manufacturer, technical educational equipments supplier, technical educational instruments, technical educational instruments exporter, technical educational instruments manufacturer, vocational training lab instruments exporter in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com