

**Product Code . ILS-H-10924**

## Ball and Ring Apparatus



### Description

---

Ring mounted on rod with a wooden handle.




The ball passes through the ring when cold but will not pass through after being heated.

Ball And Ring Apparatus is used for demonstrating thermal expansion, comprising a captive brass ball of 19mm secured to a mounted brass ring by a chain.

Contact jLab for your Educational School Science Lab Equipments. We are best school science equipments suppliers, school science kits, school science lab equipment exporter, school science lab equipment manufacturer, school science lab equipment supplier, school science lab equipment suppliers in india.

---

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com  
 [www.indianlabsuppliers.com](http://www.indianlabsuppliers.com)