

## Description

---

Grease should be applied on the flanges to make it air tight.

Plain Dessicators are physically the same as Vaccum Desiccator.

Only difference is the absence of stop cock.


**Material:** Polypropylene/ Polycarbonate

Code No.	Size (mm)	Flange O.D. (mm)	I.D. (mm)
ILS1	150	171	150
ILS2	200	232	199
ILS3	250	276	240
ILS4	300	340	300

Contact jLab for your Educational School Science Lab Equipments. We are best educational laboratory instrument suppliers, educational laboratory instruments exporters, educational laboratory microscopes, educational laboratory suppliers, educational physics lab equipments, educational science equipments in ambala, india.

---

Indian Lab Suppliers,

Direct Contact Details   
+91-8569909696 