



Product Code . ILS-MDE-11536

Disposable Wooden Tongue Depressor

Description

Instructions for use:

Item used for depressing the tongue during examination of the patient's mouth and the throat.

Safety process:

This item is for single use only.

Break the depressor tongue immediately after use to prevent re-use (risk of cross contamination).

Specifications:

Wooden tongue depressor, with rounded extremities.

Single use.

Thickness: Approximately 1.6mm

Non sterile.




Size: approx. Size: 17 x 149mm

Conditions for stock:

Store in a clean and safe environment: avoid dust, and other environmental risk of damage.

Contact jLab for your Educational School Science Lab Equipments. We are best scientific educational lab equipments, scientific equipments manufacturers, scientific lab equipment manufacturers, scientific lab plasticware, scientific laboratory equipment manufacturer, scientific laboratory equipment supplier in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com