



**Product Code . ILS-MU-11545**

## Plastic Wash Bottle 250ml

### Description

---

#### Instruction for use:

Rinse the bottle thoroughly before filling it again.

Between two dilutions, the bottle must be cleaned and disinfected with a 0.5% detergent-disinfectant solution for medical devices (i.e. 25 ml = 1 stroke of dosing pump for 5 litres water).

Write the content, the dilution date and the expiry date on the bottle.

Make sure the seal is watertight.

#### Specifications:

**Material:** Plastic.

**Capacity:** 250 ml.

**Container:** Wash bottle, round, mouth with smaller diameter than body, with screw closure and riser tube. Bottle assembled.




Contact jLab for your Educational School Science Lab Equipments. We are best science educational equipments, science educational instruments, science lab equipments manufacturer, science

---

laboratory instruments exporter, scientific educational equipments, scientific educational lab equipment in india,

---

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com  
 [www.indianlabsuppliers.com](http://www.indianlabsuppliers.com)