


Product Code . ILS-K-11677

Liquid Accelerometer



Description

The water creates a line graph in response to the acceleration.

The wooden base fits nicely in our Dynamics Carts for studying linear motion.




Add colored water to this transparent cell and align it with linear or centripetal acceleration.

You can also hang the cell from strings and study acceleration in rotating or oscillating systems.

Food coloring and instructions included. 10.5"x5.5"

Contact jLab for your Educational School Science Lab Equipments. We are best science educational equipments, science educational instruments, science lab equipments manufacturer, science laboratory instruments exporter, scientific educational equipments, scientific educational lab equipment in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com