

Product Code . ILS-FAM-11690




Lever Lab

Description

Indian Lab Suppliers Attach the center pivot to the Workshop Stand and use Hooked Masses to create equal torques on each end and establish equilibrium. Or move the pivot near one end and find the effective torque caused by the weight of the lever itself. A complete study of balance and equilibrium. 65cm long.

Contact jLab for your Educational School Science Lab Equipments. We are best school science equipments suppliers, school science kits, school science lab equipment exporter, school science lab equipment manufacturer, school science lab equipment supplier, school science lab equipment suppliers in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com