

Product Code . ILS-FAM-11691

Loop the Loop



Description

Place the ball at the top of the track. Release the ball.

The ball travels down the track in a loop without leaving the track.

This is fun and easy and sure to get the students engaged.

Show that the ball will not have enough energy to complete the loop if its initial height is too low.

7.5" Diameter.

20" Tall.

16.25"x3.875" Base.

1" neoprene rubber ball included.

Contact jLab for your Educational School Science Lab Equipments. We are best school science lab equipment suppliers in india, school science lab equipments exporter, school science lab equipments manufacturers, school science lab equipments manufacturer, school science lab equipments supplier, school science lab equipments suppliers in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com