

Product Code . ILS-SAR-11741

Singing Rods



Description

A must for explorations of sound and waves, these rods are an easy way to demonstrate longitudinal waves as opposed to transverse waves.

They're ideal for teaching about nodes and anti-nodes in standing waves.

Our Singing Rods make a great addition to your lessons on sound and waves.

These magical rods will "sing" for you if you know the secret. Simply hold a rod between your thumb and forefinger at one of the three marked nodes.




Then rub that hand along the rod from middle to end.

And they're a great follow-up activity to see if students really understand sound and waves.

Place a little of the included special rosin on your other hand.

You get two rods (75cm and 50cm), so you can demonstrate the relationship of pitch to length.

Contact jLab for your Educational School Science Lab Equipments. We are best educational lab equipments manufacturers, educational lab instruments manufacturer, educational laboratory glassware exporter, educational laboratory microscopes, educational laboratory suppliers, educational scientific lab equipments in india.

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com