

**Product Code . ILS-SAR-11745**

## Resonance Bowl



### Description

---

With practice, you should be able to produce a standing wave, with spouts of water shooting from the four anti-nodes spaced along the entire rim of the bowl.

This seemingly simple bowl takes just the right touch to turn into a fun and engaging demonstration of the behavior of waves and their interactions.

Every rub of the handle changes the frequency of the waves in the bowl.

It is quite the memorable demonstration.

The water should spring from the mouths of the four fish engraved at the bottom of the bowl and can be identified by the four streams of water engraved on the sides of the bowl.

Contact jLab for your Educational School Science Lab Equipments. We are best maths lab equipment manufacturers, physics lab equipments manufacturers, school educational equipments, school lab equipment manufacturer, school lab equipments manufacturers scientific educational lab equipment in india.

---

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  [sales@indianlabsuppliers.com](mailto:sales@indianlabsuppliers.com)

