

Product Code . ILS-MAM-11792

Paramagnetic Putty



Description

Stretch out a small piece and use the included neodymium magnet to make it jump to the magnet.

By itself, this magnetic putty acts like just about any other putty.

You can bounce it or tear it.

It even oozes and stretches.

You can also temporarily magnetize the putty by exposing a small piece of it to a magnet and use it to attract smaller metal objects.

The stronger the magnet, the more cool stuff you can do.

Putty weighs 91g (3.21oz).

But this putty is infused with micron-sized, iron-based particles that 'come alive' when exposed to a strong magnetic field.

Place the magnet on top of the putty, then watch it absorb the magnet like "The Blob!"

Contact jLab for your Educational School Science Lab Equipments. We are best scientific laboratory equipments manufacturer, scientific laboratory equipments supplier, technical educational equipment exporter, technical educational equipment manufacturer, technical educational equipment supplier, technical educational equipments exporter in india.

Indian Lab Suppliers,

Direct Contact Details  +91-8569909696  sales@indianlabsuppliers.com
 www.indianlabsuppliers.com